
Frequently Asked Questions

Whole Exome Sequencing and Analysis

Q1. What is Whole Exome Sequencing ?

A1. Whole Exome Sequencing (WES) is an efficient strategy to selectively sequence the coding regions (exons) of a genome, typically human, to discover rare or common variants associated with a disorder or phenotype [1, 2]. By focusing sequence production on exons, which represents ~2.5% of the human genome, many more individuals can be examined at significantly reduced cost and time compared to sequencing their entire genomes. The most common methods rely on hybridization by oligonucleotide probes to 'capture' targeted DNA fragments, thereby enriching for exonic sequences. Targeted exonic sequences include well-established annotated coding and non-coding exons. Regions not within close proximity, on the order of 125-bases, of the targeted regions are not sequenced. Therefore, variants within introns, promoters or inter-genic regions are generally not detected.

Note, DNA samples derived from living humans must be consented for WES before acceptance at NISC for sequencing.

For more detailed information about WES, please watch this video tutorial:

<http://www.genome.gov/27545880>

Q2. How is WES performed at NISC ?

A2. NISC currently employs a solution-based probe hybridization protocol to capture (enrich for) exonic sequences from the DNA sample. The nucleic acid probes are components of a standard commercially available kit optimally designed and prepared for this purpose by the vendor [3]. The current whole-exome capture kit used at NISC is the Nimblegen SeqCap EZ Exome +UTR Library. In brief, the DNA is sheared mechanically, targeted fragments captured by probe hybridization, and then amplified before sequencing on an Illumina HiSeq 2500 instrument. NISC continually evaluates improvements in these technologies, and implements those that represent reduction in cost and time or increase in exon coverage.

Q3. What material should I send for WES ?

A3. We need a minimum of 1.5 µg of highly-purified genomic DNA (5 µg preferred) in a volume of 120 µl or less for WES. Samples should be submitted in 1.5-1.7 ml microfuge tubes (example: VWR cat. no.89000-028) or 2 ml screw cap tubes (example: Sarstedt cat. no. 72.694.007). Please DO NOT send samples in 0.5 or 0.2 ml tubes. To ensure that each sample is uniformly pure and free of infectious agents, we require that all DNAs be phenol:chloroform extracted before submission. A simple protocol is available from NISC. Ref: www.nisc.nih.gov/docs/gDNA_submission_exome_cc.pdf

Frequently Asked Questions

Q4. How should the DNA be qualified ?

A4. The investigator must submit an image of an analytical agarose gel as evidence the DNA is of good integrity, i.e., not just a 'fuzzy blob' of low molecular weight. We highly recommend Qubit for quantitation of the DNA sample, since it uses a double-strand DNA-specific method. UV absorption methods, e.g., using a NanoDrop spectrophotometer, can drastically overestimate the concentration of DNA due to RNA and small molecule contamination.

Q5. How long are the reads for WES analyses ?

A5. Typically, NISC generates read lengths of 125 bases on a HiSeq. Paired-end reads generate a total of 250 bases of sequence (125b from each end) from each fragment in the library.

Q6. How many reads are required for WES analyses ?

A6. Currently, we target 38 million paired-end 125 base reads. This amount of reads is usually sufficient to meet our minimum coverage requirement to call genotypes with an MPG score of 10 for at least 85% of the targeted bases [4, 5].

Q7. How are variants called in WES analyses ?

A7. Sequence reads produced for a sample are aligned to the human reference sequence and the results stored in BAM format. A custom analysis program, MPG (Most Probable Genotype), processes this information using a probabilistic Bayesian algorithm, calling genotypes at all reference positions at which there are high quality bases from the aligned sequence reads [4]. The likelihood of each possible genotype from the observed sequence data is calculated and given an MPG score, where MPG 10 is considered accurate. These genotype calls have been compared against Illumina Human 1M-Quad genotype chips, and genotypes with a MPG score of 10 or greater show >99.89% concordance with SNP chip data [4].

Q8. What data are returned by NISC ?

A8. All variants, genotypes, and annotations are delivered to the investigator in tab-delimited format compatible with VarSifter [6], a java-based genotype viewer, available from NISC. The file can also be imported to Excel. The VarSifter file contains all discovered variants with genotypes of all samples sequenced, as well as gene locations (5' UTR, 3' UTR, coding-synonymous, nonsynonymous, or stop, splice site, or intron).

Frequently Asked Questions

References:

1. Biesecker L (2010) “Exome sequencing makes medical genomics a reality.” *Nature Gen.* 42, 13-14.
2. Illumina (2013) “An Introduction to Next-Generation Sequencing Technology.” www.illumina.com/documents/products/Illumina_Sequencing_Introduction.pdf
3. NimbleGen (2013) “SeqCap EZ Library.” www.nimblegen.com/products/lit/05227887001_SeqCapBroch_Jan2013.pdf
4. Teer, JK *et al.* (2010) “Systematic comparison of three genomic enrichment methods for massively parallel DNA sequencing.” *Genome Res.* **20**: 1420-1431
5. Sims, D., *et al.* (2014) “Sequencing depth and coverage: key considerations in genomic analyses.” *Nature Rev. Genetics* **15**: 121-132.
6. Teer, JK, *et al* (2012) “VarSifter: Visualizing and analyzing exome-scale sequence variation data on a desktop computer.” *Bioinformatics* **28**: 599-600.